
DentaSonic
Swiss Proximal Restorative Clean-up Tools

Products

Instructions
for use
and safety
recommendation

GB D I F E P

Product description

- Double sided diamond-coated oscillating metal strip
- Strip height: 3.6 mm
- Length of the diamond-coated part: 13 mm
- Sterilizable and reusable

Use in combination with:

- DentaSonic, W&H, KAVO, NSK, SIRONA or other oscillating HP with or whiteout ligh, and with or whiteout water irrigation.
- Manual use with handle

Indications for Use

- DiaStrip is indicated for the removal of set cement, to break contacts, for Interproximal Reduction (IPR) procedures, and for polishing and finishing.
- DiaStrip (single-sided, double-sided and resin) is indicated for both anterior and posterior teeth.
- Pass the strip back and forth gently and buccolingually only.

Contraindications

DiaStrip is contraindicated if the Instructions for Use (IFU) is not strictly followed or with patients with poor enamel structure (hypoplasia), high caries index, poor oral hygiene or severe tooth crowding.

Warnings

- Avoid bending and levering: use of excessive force (strong levering) could lead to the breakage of the instrument.
- Avoid jamming or levering actions when oscillating, as this increases the risk of instrument breakage.
- Avoid temperatures above 180°C which may affect the durability of the instrument.
- Avoid applying forces greater than 3N as this could cause damage to the instrument.
- To prevent patient from ingesting or inhaling the DiaStrip, prior to using DiaStrip, thread a piece of dental floss through the arc of the DiaStrip. Secure by tying a slipknot. Wrap the floss around one of your other fingers to secure the DiaStrip.
- The consumption of the instrument lifetime until its worn out, replacement and disposal is an exclusive decision of the physician according to his operative needs.
- To ensure traceability during the instrument lifetime we recommend to preserve the packaging.

Principles of safety

- Contact pressure: bend the strip 1-1.5 mm to achieve the ideal 1-2.5N operating force.
- Do not stop the motion of the file while the same is in the interdental space
- In case of overcrowding open the interdental space using a wedge in the cervical area
- Do not bend, warp or leverage the instrument
- Not observing these precautions leads to a high risk in terms of safety and can lead to the possible breakage or damage of the instrument

Precautions

- Sterilize DiaStrip prior to patient use. See details below.
- Read the Instructions for Use (IFU) completely before using. Use of DiaStrip inconsistent with IFU is at the sole discretion and sole responsibility of the user.
- Do not use any DiaStrip that is damaged, bent, broken, or show signs of wear.
- For maintaining anatomical contour in the posterior and to refine the natural anatomy of the anterior, finish with the curved version.
- Use a gap gauge frequently during the procedure when gap width accuracy is critical.

Adverse Reactions

- None known. However, DiaStrip is made from stainless steel and contains nickel. Nickel may cause allergic reactions in some patients.

No Reactions for allergic patients

- Resin Diamond DiaStrip is made from resin and **NOT** contains nickel. Resin will **NOT** cause allergic reactions to patients.

Instructions for Cleaning, Sterilization and Disposal

Warnings: DO NOT use cold sterile or dry heat sterilization processes.

Limitations on reprocessing: DiaStrip can be sterilized up to 10 times. Resin DiaStrip can be sterilized up to 3-5 times.

INSTRUCTIONS

Point of use: Sterilize before first use and between patients.

Preparation for decontamination: It is best to wipe immediately after use to avoid contaminants drying, thereby making the cleaning process more time-consuming and difficult. Immerse or soak DiaStrip in a hospital grade enzymatic detergent such as enzymatic solution for 15 minutes.

Cleaning: Manual Clean thoroughly - Scrub all external surfaces with a soft bristle brush until all visible soil has been removed. Use a small diameter brush cleaner to clean cannulation holes. Inspect for visible soil on exposed surfaces. Ultrasonically clean using a solution that is recommended for dental instruments such as an enzymatic solution for 5 – 10 minutes. Rinse DiaStrip under copious amount of warm water for 2-3 minutes and visually inspect to confirm removal of debris.

Do not use highly aggressive chemical products (e.g. hydrochloric acid, hydrogen peroxide), as they could corrode the instruments. **These substances can also impair the technical properties of the plastics, i.e. change their hardness and durability.**

Cleaning: Automated Not recommended

Disinfection: Not recommended

Drying: Allow to air dry in clean area. Blow lumens with clean air using filtered air source or syringe.

Maintenance, inspection and testing: Visually inspect for fractured or worn abrasive strip that could prevent proper function. Discard in Sharps container.

Packaging: Bag in an FDA cleared sterilization pouch, Sterilization pouches are recommended.

Sterilization: EU and US Autoclave Parameters

Steam autoclave- Gravity (Wrapped)

Minimum Exposure Temperature: 121°C (250°F)

Minimum Exposure Time: 30 minutes

Minimum Dry: 30 minutes

EU Autoclave Parameters

Steam autoclave- Gravity (Unwrapped)

Minimum Exposure Temperature: 134°C

Minimum Exposure Time: 3 minutes

Minimum Dry: None

US Autoclave Parameters

Steam autoclave- Pre-Vacuum (Wrapped)

Minimum Exposure Temperature: 135°C (275°F)

Minimum Exposure Time: 3 minutes

Minimum Dry: 30 minutes

Storage: DiaStrip should be stored in the sterilization pouch until ready for use with a maximum storage period of 6 months before product needs to be re-sterilized. Unwrapped DiaStrip should be used immediately after sterilization.

Recommendations for use

- Apply DentaSonic DiaStrips in interproximal areas according to the IPR enamel reduction protocol for orthodontic treatments.
- A pre-opening of tight contact points using DentaSonic DiaStrips Saw is suggested, Ref. DSS-1
- In case of very tight contact points separate the teeth by placing a wedge cervically. This makes treatment more comfortable both for operator and patient.
- Activate water spray during the entire treatment (minimum 50 mL/min).
- Regulate water flow in case of water excess or shortage.
- Insert the strip in a cervical direction while the contra-angle is moving. This prevents the instruments from jamming.
- Once the tool is inserted into the hand-piece head be sure it is fully and firmly fixed
- Switch on the hand-piece by turning on the micro-motor
- Set and bring the instrument to the recommended speed before applying it on the tooth
- With the file in motion enter the inter-proximal space to be treated with a cervical motion
- Proceed with the treatment paying attention to the speed and the pressure applied

Risk Warnings

- Avoid bending and levering: use of excessive force (strong levering effect with the contra-angle) could lead to the breakage of the instrument.
- Avoid jamming or levering actions when oscillating, as this increases the risk of instrument breakage.
- Do not exceed the 10'000 oscillations per minute to avoid breakage of the instrument.
- Avoid temperatures above 180°C which may affect the durability of the instrument.
- **Avoid applying forces greater than 3N** as this could cause damage to the instrument.
- The consumption of the instrument lifetime until its worn out, replacement and disposal is an exclusive decision of the physician according to his operative needs.
- To ensure traceability during the instrument lifetime we recommend to preserve the packaging.

Instructions for Use

Cement Removal and Break Contacts

Use DiaStrip Saw (Blue) serrated strips to remove set cement and break contacts. • Insert the serrated blade between the teeth and with a gentle "sawing motion", to remove the cement or to break contacts.

Finishing and Polishing

Use DiaStrip Single-Sided Abrasive Strips for both anterior and posterior teeth. When following the contour of teeth is critical, use the Resin Abrasive Strips as they are better suited to this purpose. • Select the desired grit metal strip (White – SuperFine, Yellow – Fine). • Select the desired grit Resin Strips (Blue – UltraFine, Yellow – Superfine, Green – Fine, White – Medium or Gray – Coarse). • Hold the DiaStrip or the Resin DiaStrip between your thumb and index finger. • In an angled position, begin a gentle "sawing motion" until light resistance is felt. • If needed, continue to other grits until the proper finish/polish is achieved.

Interproximal Reduction (IPR)

Use Single-sided and Double-Sided DiaStrip in the sequence shown below for IPR procedures. When following the contour of teeth is critical, use the Resin DiaStrip Strips as they are better suited to this purpose. For creating IPR spaces from 0.08mm to 0.2mm, use the Single-Sided DiaStrip Abrasives; for IPR spaces greater than 0.2mm, use the Double-Sided DiaStrip Abrasives. • Regardless of the strip type chosen (single, double-sided, straight metal or Resin) as described above, begin with the superfine grit (White). • Hold between the thumb and index finger. In an angled position, begin a gentle "sawing motion" until very light resistance is felt. • It is best to check the gap with a verified gap measuring gauge before moving to the next color (grit) in the sequence. Once the desired gap is achieved there is no need to continue to the next larger grit size. • If additional space is needed, switch to the fine grit (Yellow) and use in the same manner as described above. Check gap width. • If additional space is needed, switch to medium grit (Red). Check gap width. • If additional space is needed, switch to the coarse grit (Blue). Check gap width. Continue until the desired gap is achieved. • If needed, continue to lesser grits until the proper finish/polish is achieved.

IPR Sequence Guide For Orthodontic Procedures

In order to create precise interproximal spaces for Orthodontic procedures, use DiaStrip color coded strips according to the sequence outlined below. The reduction process (IPR) should begin with the white single-sided grit (0.08 mm), and continue in sequence until the desired space is achieved. Check the reduction with an approved gap measurement gauge before moving to the next color (grit) in the sequence. If needed, continue to lesser grits until the proper finish/polish is achieved.

Note: Use all versions of DiaStrip (single-sided – ESM1, double-sided – ESM, and Resin – ES) in any area of the mouth. For maintaining anatomical contour in the posterior and to refine the natural anatomy of the anterior, finish by using the Resin Diamond coated DiaStrip version.

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