



## NATURAL NITI

MBT\*\*\* SQUARE AND ROTH COMPATIBLE




## PROGRESSIVE NITI

PROGRESSIVE XL IS DAMON® COMPATIBLE




## NATURAL STAINLESS STEEL

MBT\*\*\* SQUARE AND ROTH COMPATIBLE




## PROGRESSIVE STAINLESS STEEL

PROGRESSIVE XL IS DAMON® COMPATIBLE




Archforms pictured are actual size

**BROAD NITI**



**YOUTH**

MBT\*\*\* OVOID COMPATIBLE.



**Stainless Steel**

**Permanent Electro-Etched Midlines**


Upper	3 marks
Lower	1 mark
Progressive XL	2 marks
Natural Mid-size	2 marks

**NiTi**

**Ink Midlines**

Upper – Black
Lower – Red
Progressive XL – Blue
Natural Midsize – Yellow

**LINGUAL UNIVERSAL**



\* Highland Metals is in no way affiliated with or endorsed by Ormco Corporation or Dr. Dwight Damon. (Damon and Damon SL are trade names for SDS, Ormco, and Sybron Dental Specialties)

\*\* Highland Metals prescriptions are not implied to be an exact version of any other system. Nor do we claim any endorsement by the doctor.

\*\*\* MBT is a trademark of 3M Unitek Corporation.

Visit [www.highlandmetals.com](http://www.highlandmetals.com)  
for more available arch forms.