

Titanol® Low Force

Individually packed




Description	Straight-Arch-Form		Euro-Smile-Form		Profile	ø mm	ø inch	Content
	Order No. Maxillary	Order No. Mandibular	Order No. Maxillary	Order No. Mandibular				
Titanol® Low Force with 30 % less restoring force than Titanol Superelastic up to 30 % less friction	280-0840	280-0940	-	-		0,40	.016"	10
	280-1440	280-1540	281-1440	281-1540		0,41 x 0,41	.016" x .016"	10
	280-2040	280-2140	281-2040	281-2140		0,41 x 0,56	.016" x .022"	10
	280-2046	280-2146	-	-		0,46 x 0,64	.018" x .025"	10

Titanol® Triple Force


Straight-Arch-Form

Individually packed



The 3 force zones

With the three force segments for a biologically favorable force on the teeth. Light for anteriors, medium for bicuspid and heavy for the posterior segment.



Description	Order No. Maxillary	Order No. Mandibular	Profile	ø mm	ø inch	Content
Triple Force arch wires in Straight-Arch form	276-1440	276-1540		0,41 x 0,41	.016" x .016"	10
	276-2041	276-2141		0,41 x 0,56	.016" x .022"	10

* up to 30 % less friction